


PENILE IMPLANT TACTRA (Boston Scientific)

The Tactra penile implant comprises two semi-rigid cylinders inserted into the phallus.



Source :

<https://www.bostonscientific.com/content/dam/bostonscientific/uro-wh/general/ams/hcp-resource-center-update/Tactra-13mm-Straight-Beauty.jpg>



<https://www.bostonscientific.com/content/dam/bostonscientific/uro-wh/general/ams/hcp-resource-center-update/Tactra-illustration.jpg>

OPERATING INSTRUCTIONS

Activation:

Holding the phallus with your hand, raise the phallus to place it straight to rigidify the implant.

Deactivation:

Press the phallus downwards to fold the implant.

For more information on the product:

Consult the following video: <https://youtu.be/8fIHj6ARgD0>

Complete product information document by the manufacturer:

<https://www.bostonscientific.com/content/dam/bostonscientific/uro-mh/resources/erectile-dysfunction/tactra/MH-643303-AA-Tactra-Web-Prescriptive-Information-FINAL.pdf>